

version 1.0



# 2021-2022 BROCHURE

 **major  
science**

*Innovative Life Sciences Tools*



*Incubator*



# Incubator **MS Mini Incubator**

Major Science's Mini Incubator is designed for personal use and small laboratories, saving much of space. The unit features a broad temperature range to meet a variety of microbiology or hematology applications.



- Ideal for microbiology or hematology applications
- Corrosion resistant metal chamber
- Door with large viewing area
- Backlit colored touched panel
- One stainless steel shelf is included

Cat. No.	MO-MINI
Display	LCD
Temperature Range	Ambient +5° to 70°C
Capacity	17L
Exterior Dimension (WxLxH)	310x306x380mm
Interior Dimension (WxLxH)	261x255x255mm
Weight	Approx. 13kg
Power	100-240V~, 50/60 Hz, 2A
Material	Metal
PC Connection	USB

MO-SH2524  
Stainless steel shelf for MS Mini Incubator



## MS Oven / MS Hybridization Oven


	MO-A01	MO-AOR
Cat. No.		
Display	Touch screen & graphical interface, 3.5" 64 K color TFT display	
Rotisserie / Speed	N/A	Yes (optional) 5-100rpm
Shaker motion	N/A	Orbital clockwise / counterclockwise
Shaker speed	N/A	0-200rpm
Temperature control range	Ambient +5°C to 85°C	
Temperature uniformity	± 0.2°C at 37°C	
Temperature accuracy	± 0.2°C at 37°C	
Inner chamber dimensions	(WxLxH) 340x225x260 (mm)	
Data logging	RS-232	

\*To perform rocking motion on optional accessories, MO-SEESAW is required.

- User temperature calibration available
- Timer with alarm function
- Safety door switch device
- Large color touch screen and graphical control interface for easy access and operation
- Data logging software for monitoring and recording purposes (optional)



FC CE



• Various accessories available

\*Note each accessory do not apply to all the models of MS oven. Please visit our website at [www.majorsci.com](http://www.majorsci.com) for more details.

\*All technical specifications are subject to change without notice.

\*Please visit our website at [www.majorsci.com](http://www.majorsci.com) for more product selection and detailed information.





#### US Office

19959 Sea Gull Way  
Saratoga, CA 95070  
U.S.A.  
T/ +1-408-366-9866  
F/ +1-408-446-1107

#### Taiwan Office

**Head Office:**  
No. 156, Sec. 1, Guoji Rd., Taoyuan Dist.,  
Taoyuan City 33061, Taiwan  
T/+886-3-3762878  
F/+886-3-3761310

#### Factory:

No.19, Ln. 207, Huakang St. Bade Dist.,  
Taoyuan City 33464 Taiwan  
T/ +886-3-3623319  
F/ +886-3-3623133

#### Shanghai Office

Room 612, International business exhibition center,  
9300 Hunan Road, Pudong, Shanghai, China  
National toll-free No.:400-823-9177  
T/ +86-21-50795277  
F/ +86-21-50795277

#### India Office

D. No. 12-13-99, Satguru Apts. Extn.  
Street. No. 3, Lane No.1, Tarnaka  
Secunderabad – 500 017. India  
T/ +91-40-27001515  
T/ +91-40-27001586



[www.majorsci.com](http://www.majorsci.com)

[info@majorsci.com](mailto:info@majorsci.com)