

version 1.0



2021-2022 BROCHURE

 major
science

Innovative Life Sciences Tools



Shaker





- Heavy loading capacity (up to 15kg)
- Up to 200rpm for extreme performance
- Max. of 9999 minutes timer with alarm warning
- Interchangeable / stacking platforms and accessories for different applications
- 30x30cm shaking platform for standard laboratory practice


50~2000ml flask holder
SI-200-08~13

Platform with non-slip
rubber mat
MS-P3030

30x30cm dimpled mat
MS-DIMPLED-30

20x20cm sticky pad
PPL-04-SI-SI-200-3

Universal adjustable platform
MS-UP-30



*Anti-moisture shakers (MS-NOR-3001/MS-NRK-3001/MS-NRC-3001) are also available for cold room and CO₂ incubation operation. Please contact your regional manager for detail information.



Cat. No.	MS-NOR-30	MS-NOR-3001	MS-NRK-30	MS-NRK-3001	MS-NRC-30	MS-NRC-3001
Shaking motion	Orbital action in one direction or two directions		Rocking Tilt angle:12°		Reciprocal Stroke length:20mm	
Speed/Inc.	0 ~ 200rpm / 1rpm		5 ~ 100rpm / 1rpm			
Timer	9999 (min) with alarm / Continuous					
Platform dimension (cm) (WxL)	30x30					
Loading capacity	10kg		15kg		8kg	
Stacking platform	Yes					
Anti-moistured	N/A	Yes	N/A	Yes	N/A	Yes
Platform capability of flask holders	50ml	13			13	
	125ml	12			12	
	250ml	9			9	
	500ml	5			5	
	1000ml	4			4	
	2000ml	1			2	

◆ Customized specification



US Office

19959 Sea Gull Way
Saratoga, CA 95070
U.S.A.
T/ +1-408-366-9866
F/ +1-408-446-1107

Taiwan Office

Head Office:

No. 156, Sec. 1, Guoji Rd., Taoyuan Dist.,
Taoyuan City 33061, Taiwan
T/+886-3-3762878
F/+886-3-3761310

Factory:

No.19, Ln. 207, Huakang St. Bade Dist.,
Taoyuan City 33464 Taiwan
T/ +886-3-3623319
F/ +886-3-3623133

Shanghai Office

Room 612, International business exhibition center,
9300 Hunan Road, Pudong, Shanghai, China
National toll-free No.:400-823-9177
T/ +86-21-50795277
F/ +86-21-50795277

India Office

D. No. 12-13-99, Satguru Apts. Extn.
Street. No. 3, Lane No.1, Tarnaka
Secunderabad – 500 017. India
T/ +91-40-27001515
T/ +91-40-27001586



www.majorsci.com

info@majorsci.com